

Donastar Commercial Water Filtration

Contact us for more
information

(833) 366-2782

support@donastar.com

CTO = Chlorine, Taste, and Odors
P1 = Particulate Class 1 (particles 0.5 to 1 micron in size)
TSS = Total Suspended Solids
CLM = Chloramine

Filter Series ⁽¹⁾	Reductions / Description	Pore Size (micron)	Length (in)	Pressure (psi)	Temp. (°F)	CTO (GPM)	CTO Capacity (kGal)	CLM (GPM)	CLM Capacity (kGal)
14-Series									
70-14090-x41			3x14.75			2.0	9.0	–	–
70-14150-x61	CTO, P1, cysts; bacteria static	1	3x17	30-100	35-100	2.5	15.0	–	–
70-14220-x61			3x17			2.0	22.0	–	–
15-Series									
70-15125-x41	Scale, CTO, P1, cysts, bacteria static	1	3x14.75	30-100	35-100	2.0	12.5	–	–
70-15220-x61			3x17			2.0	22.0	–	–
30-Series									
70-30000-E61 ⁽²⁾	TSS, 5nm+ particulates	0.005 ⁽²⁾	3x17	30-75	40-100	1.5	9.0	–	–
34-Series									
70-34095-x61 ⁽³⁾	Soluble lead, CTO, P1, cysts; bacteria static	1	3x17	30-100	35-100	2.5	9.5	–	–
35-Series⁽³⁾									
70-35220-x61	Scale	N/A	3x17	30-100	35-100	2.0	22.0	–	–
70-35350-x61						2.0	35.0	–	–
70-35550-x61						2.0	55.0	–	–
36-Series									
70-36000-x61	Sediment	5	3x17	10-100	40-100	–	–	–	–
41-Series									
70-41065-x41	Chloramine, CTO, scale	1	3x14.75	30-100	35-100	1.38	9.0	1.38	6.5
70-41075-x61			3x17			1.67	12.0	1.67	7.5
50-Series									
70-50000-A61	Air and gas contaminants	5	3x17	10-80	35-100	–	–	–	–
Drop-Ins, Pre-Filters, etc.									
70-35100-S10	Drop-In, Scale	N/A	2.5x10						
SDC-25-2005	Pre-Filter, Sediment	5	2.5x20						
70-38001-S10	UV Twin System Carbon	5	–						

Manifolds

70-00000-AHF1	High Flow Single-Stage Manifold	70-00000-AHF4	High Flow Quad-Stage Manifold
70-00000-AHF2	High Flow Dual-Stage Manifold	70-00000-AST1	Smart Tap Single-Stage Manifold
70-00000-AHF3	High Flow Triple-Stage Manifold		

(1) "x" in product number indicates availability in the following keys: A = Proprietary and E = Everpure; C = Cuno key also available with minimum order quantity of 1,000 units.

(2) 30-series filters are currently available only in our E key; A and C keys may be available upon request; 70-30000 is not an FDA-cleared device, nor intended for use in infection control.

(3) 34-series and 35-series filters are currently available only in our A and E keys; C key may be available upon request.

DISCLAIMER: Everpure® is a registered trademark of Pentair PLC and Cuno™ is a trademark of 3M Company, each of which is not affiliated or otherwise associated with Nephros, Inc. Nephros, Inc. acknowledges and confirms that it has not been granted any license or authority from Pentair PLC or 3M Company to use the marks, and further confirms that any and all use of the marks on this site is for identification and descriptive purposes only and such use is not intended to, does not, and shall not be deemed to imply any association or affiliation between Nephros, Inc. and Pentair PLC or 3M Company of any nature whatsoever.

Manufactured in the USA by Nephros, Inc., 380 Lackawanna Place, South Orange NJ 07079. Donastar Enterprises and Nephros are registered trademarks. All rights reserved.

DONASTAR.COM

DONASTAR ENTERPRISES LLC:

421C Street
Suite 1A
Washougal, WA 98671 USA
(833) 366-2782

021023